

PELCO® Prep-Eze™ Rectangular Wellplate Inserts - Instructions

Product No's. 36168, 36170 and 36172



For use with Bench and Microwave Processing Protocols

The Prep-Eze rectangular wellplate inserts are designed for use with standard polystyrene tissue culture 6, 12 and 24 flat-bottom multiwell culture plates (6-well Prod. #36168-6, 12-well Prod. #36170-12 and 24-well Prod. #36172-24). Designed for the labeling of free-floating sections, on the bench or in the microwave.

Contents

The PELCO® Prep-Eze™ Rectangular Wellplate Insert (6, 12 or 24 well) is enclosed. The insert is made from solvent resistant materials and the screens are replaceable (See Fig. 1 below).



Fig. 1. The screen in each well fits into a recessed groove as shown. The screen can be removed by using a flat/blunt object like a pencil eraser to push the screen out. A new screen is put in place by inserting it into the top of the well and gently pushing down with the eraser end of a pencil or similar blunt object.

Set Up

The PELCO® Prep-Eze™ Rectangular Wellplate Insert is ready to use out of the box. For easy sample identification the insert has the same square and curved corners as the multiwell culture plates.

Use

The wellplate insert is placed in the appropriate sized multiwell culture plate and the reagent volume determined for the staining and/or processing regiment being undertaken. Each well can then be filled with the appropriate volume required for the procedure. Specimen transfers are easy and quick.

NOTE: Do not use any reagents that will damage, dissolve or adversely interact with the multiwell culture dish being used. Cleaning of the PELCO® Prep-Eze™ Rectangular Wellplate Insert can be done with hot soapy water or other similar cleaning agent. Autoclaving should be avoided. Chemical sterilization is the recommended choice.

WARNING:

The PELCO® Prep-Eze™ Rectangular Wellplate Inserts should not be used for staining or processing at temperatures over 70°C.

36168 36170 36172 TN V1 06272003

TED PELLA, INC.
Tools for Science and Industry

P.O. Box 492477, Redding, CA 96049-2477, U.S.A.
Telephone: 530-243-2200; 800-237-3526 (U.S.A. or Canada) • FAX: 530-243-3761
Email: sales@tedpella.com • Web Site: <http://www.tedpella.com>

ISO9001/9002
REGISTERED
URS Cert. 2629A/2629B 9/00

Replacement Mesh and Accessories

Product	Product Accessory No.	Description	Qty
6-well Prod. #36168	36168-8	6-well Insert Replacement Mesh	pkg/30
	36168-6	6-well cavity Tissue Culture Plate	pkg/50
12-well Prod. #36170	36170-14	12-well Insert Replacement Mesh	pkg/36
	36170-12	12-well cavity Tissue Culture Plate	pkg/50
24-well Prod. #36172	36170-24	24-well Insert Replacement Mesh	pkg/48
	36172-24	24-well cavity Tissue Culture Plate	pkg/50

TED PELLA. INC.

Tools for Science and Industry

P.O. Box 492477, Redding, CA 96049-2477, U.S.A.
Telephone: 530-243-2200; 800-237-3526 (U.S.A. or Canada) • FAX: 530-243-3761
Email: sales@tedpella.com • Web Site: <http://www.tedpella.com>

