

NEW

PELCO BioWave® Pro+


Trusted Technology, Enhanced Efficiency

The PELCO BioWave® Pro+ introduces new features that improve efficiency and streamline functionality. It is built on trusted technology and adds enhancements for efficient laboratory processing.

- + User-friendly run screens with clickable text navigation
- + Live run-time graph and clear view of run settings in real-time
- + Ideal for multi-user environments with new multi-user profiles
- + Two USB ports for easy data transfer and storage
- + Easily upload custom protocols
- + PELCO BioWave® Pro+ RPM app for data transfer and review




The PELCO BioWave® Pro+ is the trusted technology you know with the enhanced efficiency that you need.




User-Friendly Run Screens

Intuitive functions in a simple, consistent format offer all setting controls within one screen.



Live Run-Time Graph

Quick and convenient view of the current protocol in real-time provides an active view of run progress.



Simplified Protocol Selection

Icon-driven interface allows users one-touch access to pre-loaded and customized protocols.

PELCO BioWave® Pro+

Trusted technology. Enhanced efficiency.

NEW




The PELCO BioWave® Pro+ and RPM (Report/Protocol Manager) app simplify data transfer and custom protocol upload.



Two USB Ports

The PELCO BioWave® Pro+ now has two USB ports for convenient data transfer of protocols and recorded run information.



PELCO BioWave® Pro+ RPM

The Report/Protocol Manager (RPM) app is designed to meet diagnostic laboratory regulations on reviewing recorded data for each protocol. Recorded data can be downloaded to the USB key and then transferred into the RPM app on a laboratory computer.

The RPM app also provides an improved design for protocol creation and review. The completed protocol can be uploaded to the USB key and then downloaded to the PELCO BioWave® Pro+ directly from the touchscreen.

©Ted Pella, Inc. 1-5-2017, Printed in U.S.A.

www.tedpella.com

TED PELLA, INC.
Microscopy Products for Science and Industry

800.237.3526