

# BLADES FOR MATRICES & MICROTOMES



## WECPREP™ BLADES

Single edge, stainless steel, surgically sharp. Sturdy, long and effective cutting blade. Packaged in 10 small boxes of 5 blades each in larger box.

Length: 57.4mm (2.26"), Cutting Edge Length: 55.9mm (2.20"), Width: 13.4mm (0.528"), Thickness: 0.279mm (0.011")

121-10 Wecprep™ Blades, box/50.....each



## PERSONNA® STAINLESS STEEL SURGICAL PREP BLADES

Personna® single edge, stainless steel, surgically sharp blade. Smooth, long and reliable blade. Packaged in 10 small boxes of 5 blades each in larger box.

Length: 56.8mm (2.24"), Cutting Edge Length: 54.8mm (2.16"), Width: 13.3mm (0.523"), Thickness: 0.228mm (0.009")

121-13 Personna® Stainless Steel Surgical Prep Blades, box/50 .....each



## SINGLE EDGE, LONG, CARBON STEEL BLADES

High quality, single edge, long carbon steel blades for extra cutting width. Individually protected, box of 250.

Length: 57.4mm (2.26"), Cutting Edge Length: 54.9mm (2.16"), Width: 13.3mm (0.524"), Thickness: 0.229mm (0.009")

121-15 Single Edge, Long, Carbon Steel Blades, box/250 .....each



## ULTRA-THIN BLADES

Ultra-thin, single edge blade Single edge, stainless steel blade. Ultra-thin, sharp edge blade with a cutting tip angle of 22°.

Length: 39.9mm (1.55"), Cutting Edge Length: 39.9mm (1.55"), Width: 9mm (0.35"), Thickness: 0.076mm (0.003")

121-19 Ultra-Thin Blade, Single Edge, box/100.....each



## PELCO® BRAIN MATRIX BLADE

Single edge uncoated carbon steel. This single edge thin, extra long razor blade is made from uncoated carbon steel. It is double beveled for increased strength and durability. Each blade is individually protected and lightly oiled for corrosion protection. This blade is developed specifically for large brain matrices.

Length: 120mm (4.72"), Width: 20mm (0.78"), Thickness: 0.12mm (0.005")

121-96 PELCO® Brain Matrix Blade, pkg/25..... each



## BLADE SAFE

Safely store up to 100 used blades for disposal. Discard when full.

Length: 69.85mm (2.75"), Width: 50.8mm (2"), Height: 38.1mm (1.5"), Blade Slot Length: 60.325mm (2.375")

121-25 Blade Safe .....each

# BLADES FOR MATRICES & MICROTOMES




## FEATHER® BLADES

Double edge, carbon steel blade. Sharp edges; may be broken in half lengthwise (use eye and finger protection for safety). The preferred cutting blade for the PELCO easiSlicer™.

Length: 42.8mm (1.69"), Cutting Edge Length: 36mm (1.42"), Width: 22mm (0.865"), Thickness: 0.127mm (0.0045")

121-9 Feather® Double Edge Carbon Steel Blades..... box/10



## DOUBLE EDGE, BREAKABLE STYLE BLADES

PTFE-treated, individually wrapped and bulk packed. Lightly oiled to prevent corrosion.

Length: 42.9mm (1.69"), Cutting Edge Length: 37.1mm (1.46"), Width: 21.7mm (0.855"), Thickness: 0.11mm (0.004")

121-6 Double Edge Stainless Steel Cutting blades ..... box/250




## GEM® SINGLE EDGE, CARBON STEEL BLADES

Best quality 3 facet blade, single edge, aluminum back. Individually protected, 2 boxes of 100 each.

Length: 38.9mm (1.53"), Cutting Edge Length: 38.8mm (1.527"), Width: 19.6mm (0.77"), Thickness: 0.229mm (0.009")

121-1 GEM® Carbon Steel Blades..... box/200



## SINGLE EDGE, CARBON STEEL BLADES, DEGREASED

Durable 2 facet blade, single edge, steel back. Degreased and washed to remove oil and reduce particles.

Length: 38.9mm (1.53"), Cutting Edge Length: 38.8mm (1.527"), Width: 19.6mm (0.77"), Thickness: 0.229mm (0.009")

121-32 Carbon Steel Blades, 2 Facet, Degreased, NOT Individually Wrapped..... box/115



## GEM® SINGLE EDGE, STAINLESS STEEL UNCOATED BLADES

Best quality 3 facet blade, single edge, aluminum back. Individually protected, 2 boxes of 100 each.

Length: 38.9mm (1.53"), Cutting Edge Length: 38.8mm (1.527"), Width: 19.8mm (0.78"), Thickness: 0.229mm (0.009")

121-2 GEM® Stainless Steel Uncoated Blades ..... box/200



## GEM® SINGLE EDGE, STAINLESS STEEL BLADES, DEGREASED

Same 3 facet blade as 121-2 with aluminum back. Degreased and washed to remove oil and reduce particles.

Length: 39.6mm (1.56"), Cutting Edge Length: 39.6mm (1.56"), Width: 19.6mm (0.77"), Thickness: 0.229mm (0.009")

121-31 GEM® Stainless Steel Blades, Degreased, NOT Individually Wrapped ..... box/115



## GEM® SINGLE EDGE, PTFE-COATED BLADES

Cuts through epoxies nicely. Stainless steel single edge, 3 facet, PTFE-coated stainless steel GEM® blades. Individually protected, 2 boxes of 100 each.

Length: 39.6mm (1.56"), Cutting Edge Length: 38.8mm (1.527"), Width: 19.6mm (0.77"), Thickness: 0.229mm (0.009")

121-3 GEM® Single Edge, PTFE-Coated, Stainless Steel Blades ..... box/200



## INJECTOR-STYLE, STAINLESS STEEL, PTFE-TREATED BLADES

Useful for microtome sectioning and other sectioning applications. Single edge, 20 per dispenser.

Length: 38mm (1.5"), Cutting Edge Length: 36.3mm (1.43"), Width: 8mm (0.314"), Thickness: 0.25mm (0.010")

121-4 Injector Blades, PTFE-Treated..... dispenser/20

