



CellSafe™+ Biopsy Cassette v2 Secure Biopsy Specimen Handling

Benefits

- Improved security, new locking latch
- Firmer mesh
- Blue mesh for better tissue contrast
- X-Ray transparent for breast work
- 20 times less carry over than foam pads
- Zero artifacts associated with foam pads

The CellSafe™+ Biopsy Capsule offers high specimen security designed to contain small tissue fragments and biopsies during tissue processing. The CellSafe™+ Biopsy cassette fits into a standard tissue cassette for secure processing. The CellSafe+ is designed for single use.

Interior dimensions of a regular cassette are: 30 x 26 x 4.92mm H.

Exterior dimensions of the CellSafe™ Biopsy Cassettes are: 27.5 x 25 x 4.57mm H.

IVD CE In Vitro Diagnostic Medical Device

27165	Cassettes, CellSafe™, Biopsy, White 27.5 x 25 x 4.57mm H.....	pkg/100
7165-8	Cassettes, CellSafe™, Biopsy, Blue 27.5 x 25 x 4.57mm H.....	pkg/100

NEW CellSafe™ 5 Biopsy Cassettes

- Five chamber biopsy cassette
- Small needle core tissue specimens may be processed separately
- Ideal for breast and prostate core biopsies
- X-ray transparent fine mesh
- Available in white and blue colors



The CellSafe™ is designed to contain small tissue fragments and biopsies during tissue processing. The CellSafe™ 5 was developed for the simultaneous preparation of core needle biopsies. Each well is numbered for identification. The blue cassette aids with contrast and visualization of light colored biopsy tissues.

If you are using any unusual reagents it is suggested that you carry out preliminary chemical compatibility tests before using the CellSafe™ 5 to contain tissue samples.

IVD CE *In Vitro* Diagnostic Medical Device

27165	Cassettes, CellSafe™, Biopsy, White 27.5 x 25 x 4.57mm H.....	pkg/100
7165-8	Cassettes, CellSafe™, Biopsy, Blue 27.5 x 25 x 4.57mm H.....	pkg/100