

LIGHTED EMBEDDING FORCEPS / SCISSOR



5003-1 LED Lighted Ergonomic Embedding Forceps 14cm Kit, Smootheach



5003-5 LED Lighted Embedding Forceps 15cm Kit, Smootheach



5003-6 LED Lighted Metzenbaum Scissors 14.5cm Kit..... each

HISTOLOGY FORCEPS

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS

www.tedpella.com

AESCULAP® DISSECTION FORCEPS


Finest Quality Forceps – Small/Medium with Smooth or Serrated Tips



5002-5 Aesculap® Iris Dissection Forceps, 97mm L (3.8")each




5002-7 Aesculap® Dissection Forceps, Very Delicate, 115mm L (4.5")each



5002-9 Aesculap® Adson Dissection Forceps, Curved, 120mm L (4.75")each



5002-11 Aesculap® Dissection Forceps, Smooth Jaws, 185mm L (7.25")each



5002-22 Aesculap® Dissection Forceps, Smooth Jaws, 200mm L (8")each



5002-20 Aesculap® Delicate Dissection Forceps, Stainless Steel Satin Finish, Straight, 100mm L (4")each



TED PELLA, INC.

Microscopy Products for Science and Industry

www.tedpella.com

sales@tedpella.com

800-237-3526

AESCULAP® DISSECTION FORCEPS

Medium, Serrated Tips



5002-21 Aesculap® Delicate Dissection Forceps, Stainless Steel Satin Finish, Half-Curved, 100mm L (4").....each



5002-23 Aesculap® Medium Dissection Forceps, Straight, 115mm L (4.5")each



5002-25 Aesculap® Medium Dissection Forceps, Straight, 130mm L (5.12")each



5002-27 Aesculap® Medium Dissection Forceps, Straight, 145mm L (5.71")each



5002-29 Aesculap® Medium Dissection Forceps, Straight, 160mm L (6.30")each



5002-35 Aesculap® Medium Dissection Forceps, Curved, 130mm L (5.12").....each



5002-37 Aesculap® Medium Dissection Forceps, Curved, 145mm L (5.71").....each

Your source for high quality and economy forceps & tweezers for every application.



SUTURE REMOVAL KIT

Disposable, sterile suture removal kit offers convenience and extra protection against cross infection. This sterile kit can also be used in the procurement of surgical, forensic, and pathology tissue biopsies. Sold in packs of 10.

Kit contains littauer scissors, 4" forceps and gauze sponge

13209 Suture Removal Kit, sterile, single use pkg/10



PELCO® TWEEZER TRAP

A safe convenient and portable tweezer holder in 3 different models for 3, 5 or 8 tweezers or forceps. Made from strong white acrylic with stainless steel pins. Suitable for tweezers up to 140mm length.

* Tweezers sold separately

5689 PELCO® Tweezer Trap, 3 Tweezers Capacity each
 5690 PELCO® Tweezer Trap, 5 Tweezers Capacity each
 5691 PELCO® Tweezer Trap, 8 Tweezers Capacity each



5168 Neoprene O-Rings, for Tweezer Clamping, pkg/10 each
 5169 Plastic Tweezer Tip Guards, Straight, pkg/10 each



TWEEZER REPAIR KIT

Salvage damaged tweezers with this kit. Use for straightening, sharpening and deburring dull or damaged tweezers and other small precision tools. Contains Arkansas Whetstone, Lapping Paper Strips, Neoprene O-ring Tweezer Clamps and Steel Straightening Block. Instructions included.

54528 Tweezer Repair Kit each

TWEEZER KITS FROM TED PELLA, INC.



PELCO® PRO SUPERALLOY, HIGH PRECISION TWEEZER SET

Set comes in a top quality, convenient, hard plastic wallet with foam lining.

- 5385-50SU PELCO® Pro Superalloy Tweezer Set, Styles 1, 2A, 3C, 5 and 7, in a wallet.....each
- 5367-50NM PELCO® Pro High Precision NM-SS Tweezer Set, Styles 00, 2A, 3C, 5 and 7, in wallet.....each



PELCO® DISSECTION KIT

Surgical instruments for the teaching or research laboratory. All instruments are made of durable stainless steel and have a fine finish. Supplied in a zipper case with an 8 instrument capacity.

- 13200 PELCO® General Biology Dissecting Kit pkg/10




INDUSTRIAL, ECONOMY TWEEZERS


Made from anti-acid, non-corrosive, non-magnetic stainless steel. Fine brushed finish, a great value for the price.

- 58098-NM Set of 4 Economy Tweezers, #2A, #3, #4 and #7, Stainless Steel, Anti-Acid, Non-Magnetic, in a convenient Plastic Poucheach
- 58099-NM Set of 4 Economy Tweezers, #2A, #3, #4 and #7, Stainless Steel, Anti-Acid, Non-Magnetic, in a convenient Plastic Poucheach


Your source for high quality and economy forceps & tweezers for every application.



5002-45 Standard Dressing Forceps, 254mm (10"), 3.7mm Tip, Straight, Serratedeach



5002-46 Standard Dressing Forceps, 300mm (12"), 4.5mm Tip, Straight, Serratedeach



5002-38 Semken Dressing Forceps, Fine, 127mm (5"), 1mm Tip, Straight, Serrated.....each




5002-39 Semken Dressing Forceps, Fine, 152mm (6"), 1.6mm Tip, Straight, Serrated.....each

SPLINTER FORCEPS



5002-50 Splinter Forceps, 114mm (4.5"), 1mm Tip, Straight, Smooth.....each



5002-51 Splinter Forceps, 114mm (4.5"), 1mm Tip, Straight, Smooth.....each



5002-52 Splinter Forceps, 114mm (4.5"), 1mm Tip, Straight, Serrated.....each




5002-53 Splinter Forceps, 114mm (4.5"), 1mm Tip, Curved, Serratedeach

SEE PRODUCT DETAILS AND COMPLETE SELECTIONS | www.tedpella.com/tweezers_html

HISTOLOGY FORCEPS FROM TED PELLA, INC.

DEBAKEY TISSUE FORCEPS



13221 DeBakey Tissue Forceps, 160mm (6.3”), 1.6mm Tip, Serrated.....each




13222 DeBakey Tissue Forceps, 160mm (6.3”), 1.5mm Tip, Serrated.....each




13223 DeBakey Tissue Forceps, 160mm (6.3”), 2.0mm Tip, Curved, Serratedeach


STANDARD DRESSING FORCEPS




5002-42 Standard Dressing Forceps, 114.3mm (4.5”), 2.5mm Tip, Straight, Serratedeach



5002-43 Standard Dressing Forceps, 152.4mm (6”), 3.0mm Tip, Straight, Serratedeach



5002-41 Standard Dressing Forceps, Delicate, 203mm (8”), 3.0mm Tip, Straight, Serratedeach



5002-40 Standard Dressing Forceps, 203mm (8”), 3.5mm Tip, Straight, Serratedeach

Your source for high quality and economy forceps & tweezers for every application.



13220 Tissue Forceps, 1 x 2 Teeth, 127mm L (5")each

STAINLESS STEEL MICROSCOPY FORCEPS



5441 Microscopy Forceps, 90mm L (3.54")each



5443 Microscopy Forceps, 115mm L (4.53")each



5445 Microscopy Forceps, 145mm L (5.71")each

ECONOMY DISSECTION FORCEPS

Curved or Straight with Serrated Tips



13250 Straight Dissection Forceps, Serrated Tips, 127mm L (5")each



5431 Curved Dissection Forceps, Serrated Tips, Chrome Steel, 117mm L (4.6")each



5432 Curved Dissection Forceps, Serrated Tips, Stainless Steel Polished Finish, 140mm L (5.5"), 0.4 x 0.4mm Pointeach

SEE PRODUCT DETAILS AND COMPLETE SELECTIONS | www.tedpella.com/tweezers_html