

LARGE SLIDES & COVER GLASS



LARGE GLASS MICROSCOPE SLIDES

Especially made for larger tissue sections, such as whole brain sections, to enable uninterrupted examination w/o loss of tissue. Standard thickness is 1.2mm, cut ends (not available with frosted ends). Made from optical grade soda lime glass. Sizes available are 4" x 3" up to 7" x 5". Packed as 36 slides per shrink-wrapped package. Pre-cleaned.



260230	Large Glass Slides, 4" x 3" (102 x 76mm)pkg/36
260231	Large Glass Slides, 4" x 3¼" (102 x 83mm)pkg/36
260232	Large Glass Slides, 5" x 4" (127 x 102mm)pkg/36
260233	Large Glass Slides, 6" x 4½" (152 x 114mm)pkg/36
260234	Large Glass Slides, 7" x 5" (178 x 127mm)pkg/36



LARGE PLAIN ECONOMY MICROSCOPE SLIDES

Made from clear soda lime corrosion resistant glass. Precleaned slides with a size of 3" x 2" and a 0.039" - 0.047" thickness (76.2 x 50.8mm and a 1 - 1.2mm thickness). Plain, ground edges and beveled corners, 144 slides per package.

260439	Large Plain Glass Slides, Beveled Corners, 3" x 2"	pkg/144
--------	--	---------



HISTOBOND®+ SUPER MEGA SLIDES & COVER GLASS

These slides have thin white frosted end for labeling. The unfrosted area is sized to accommodate the standard large section cut from a paraffin embedded tissue block from the large Super Mega and Super Mega Slim Base Molds (60 x 45mm).

Dimensions: 75.2 x 50.4 x 1mm (3" x 2" x 1mm)

260240	HISTOBOND® + Super Mega Slides	pkg/50
260378	Coverglass for Super Mega Slides, #1.5 Thickness	pkg/100



260238	Adhesion Super Frost + Slides, 38 x 75mm (1-1/2" x 3") pre-cleaned	1 gross
260239	Adhesion Super Frost + Slides, 51 x 75mm (2" x 3") pre-cleaned	1 gross

CORNING® LARGE GLASS MICROSCOPE SLIDES

Finest quality noncorrosive glass, with finely ground edge. Pre-cleaned and packed vertically. Packed 72 slides per box, 2 boxes per package.

No notable interference from UV absorption or reflection in the 250-350nm range (tested with 280nm).

Thickness: 0.96mm - 1.06mm (0.038" - 0.042")



26005	Corning® Glass Slides, 75 x 50mm	pkg/144
-------	--	---------

LARGE COVER GLASS: THICKNESS #2

Large cover glass specially made for large microscope slides. Thickness #2 (0.19-0.25mm). Made from optical grade borosilicate glass. Sizes available 3" x 3¼" up to 4½" x 6¼". Packed as 36 cover glasses per shrink-wrapped package.

260450	Large Cover Glass 3" x 3¼" (76 x 83mm)	pkg/36
260451	Large Cover Glass 3" x 3½" (76 x 89mm)	pkg/36
260452	Large Cover Glass 3¾" x 4½" (95 x 114mm)	pkg/36
260453	Large Cover Glass 4⅛" x 5½" (105 x 140mm)	pkg/36
260454	Large Cover Glass 4½" x 6¼" (114 x 159mm)	pkg/36

LARGE COVER GLASS: THICKNESS #1 AND #1½

Made of optical grade borosilicate glass. Made for large microscope slides 3"x 2" or larger. Thickness #1 (0.13-0.17mm). Thickness #1½ (0.16-0.19mm).

Thickness #1

260460	Cover Glass 36 x 60mm (1.42" x 2.36"), 37 pcs/oz.	pkg/1oz.
260363	Cover Glass 43 x 50mm (1.70" x 1.97"), 35 pcs/oz.	pkg/1oz.
260364	Cover Glass 48 x 60mm (1.89" x 2.36"), 26 pcs/oz.	pkg/1oz.
260365	Cover Glass 48 x 65mm (1.89" x 2.56"), 24 pcs/oz.	pkg/1oz.
260462	Cover Glass 50 x 75mm (1.96" x 2.95")	pkg/100

Thickness #1½

260461	Cover Glass 36 x 60mm (1.42" x 2.36"), 30 pcs/oz.	pkg/1oz.
--------	--	----------

SLIDE STORAGE BOXES

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS

www.tedpella.com



WOODEN MICROSCOPE SLIDE STORAGE BOXES

Durable, hardwood construction for 100 slides. Available for 1½" x 3" or 2" x 3" (38 x 75mm or 50 x 75mm) large slides. The lid interior contains a permanent inventory sheet for cataloging and identification. Hinged lid with a secure latch lock closure.

- 21079 Wooden Storage Box for 100 Slides 38 x 75mm each
 21077 Wooden Storage Box for 100 Slides 50 x 75mm each



TALL SLIDE BOX

Made of sturdy ABS material, this box is ideal for storage and transport of tall 3" x 2" (75 x 50mm) microscope slides. The grooved slots hold the slides securely in place. The base contains a cork liner with location numbers for convenient inventory. Hinged lid with a secure latch lock closure. The lid interior contains a permanent inventory sheet for cataloging and identification.

- 2195 Tall Slide Box, Holds 25 Slides (2 x 3" slides) each
 2197 Tall Slide Box, Holds 100 Slides (2 x 3" slides) each



SOLID HARDWOOD SLIDE STORAGE BOXES

- Numbered index on inside cover
- Hinged cover with two latches
- Handle on one end for easy storage and stacking access
- High quality construction and finish
- Available in 5 sizes to accommodate all large slides

Solid Hardwood Storage Box for 100 Slides

- 2186 Solid Hardwood Storage Box 7¾" x 17⅜" x 4¾" (197 x 441 x 121mm), Holds 3" x 4" Slides each
 2187 Solid Hardwood Storage Box 8¼" x 17⅜" x 5" (210 x 441 x 127mm) Holds 3¼" x 4" Slides each

Solid Hardwood Storage Box for 50 Slides

- 2188 Solid Hardwood Storage Box 6⅜" x 17⅜" x 5" (162 x 441 x 127mm) Holds 4" x 5" Slides each
 2189 Solid Hardwood Storage Box 7⅜" x 17⅜" x 5¼" (187 x 441 x 133mm) Holds 4½" x 6" Slides each
 2190 Solid Hardwood Storage Box 8¼" x 17⅜" x 5¾" (210 x 441 x 146mm) Holds 5" x 7" Slides each



BLOCSTOR3 FOR SUPER MEGA CASSETTES AND SLIDES

The BlocStor range of tissue block storage products is manufactured from recycled fiberboard and provides a cost-efficient archiving and storage solution. All boxes are designed for strength, rigidity, long term

archiving of tissue blocks, and for maximum utilization of storage space.

Every BlocStor is supplied with a pre-printed archive label for clear identification and rapid retrieval. BlocStor boxes are available assembled or flat-packed for efficient use of your space.

Dimensions: 47 L x 17.2 W x 8.2cm H (16-1/4" L x 6-3/4" W x 3-1/4" H)

- 27136 BlocStor3 Super Mega Cassette Block and 2 x 3" Slide Storage Unit, 3 rows each



STURDY BOARD STORAGE BOXES

Each hold up to 15 slides per box with 5 separate compartments that can store 3 slides each in separate four flap enclosures. Stored slides are protected by the individual enclosures and easy access and removal can be done without damage. Boxes are

secured on the corners with enameled metal strips.

- 2191 Sturdy Board Storage Box 6⅞" x 6" x 4¾" (165 x 152 x 121mm) Holds 15 Slides 4" x 5" each
 2192 Sturdy Board Storage Box 8½" x 6" x 5¾" (216 x 152 x 146mm) Holds 15 Slides 5" x 7" each

 **TED PELLA, INC.**
 Microscopy Products for Science and Industry

www.tedpella.com sales@tedpella.com 800-237-3526