

LABORATORY MARKING PENS

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS
www.tedpella.com

MARKING PENS FOR HISTOLOGY



Available colors: ●

STATMARK™ PEN

The STATMARK™ Pen is the ideal black marker for histology and cytology labs. It has a permanent ink that resists even the harshest chemicals including: formalin, ethanol, isopropanol and xylene. Writing remains clean, clear and legible throughout the range of laboratory procedures – including processing and staining. The ink dries quickly and will not run or smear. It lasts through fixatives, baths, dips, stains, bleaches and more. Best of all, the STATMARK™ Pen will not readily dry out.

22313 STATMARK™ Pen.....box/12



Available colors: ●

MOIST MARK PLUS™ PEN

Premium slide/cassette marker. Writes smooth, dark and fine on all cassettes and all slides. Resistant to solvents, light resistant and water resistant on most surfaces. Use for cytology as a “dotting” reference marker as well.

22314 Moist Mark Plus™ Pen.....box/10

GENERAL LABORATORY AND CRYOGENIC PERMANENT MARKING PENS




Available colors: ● ●

SCIENCE-MARKER™

Alcohol and Water Resistant Cryogenic Marker – Permanent Lab Cryo-Marker™ with a Fine Tip (0.8 mm) Alcohol-Resistant on most surfaces. Preliminary testing is required. Not suitable for exposure to Acid/Ethanol mixes. Six markers per pack.

Available Colors: Black and Green

22315-0 Science-Marker™ Black.....pkg/6
 22315-1 Science-Marker™ Green.....pkg/6



Available colors: ● ● ● ● ● ● ● ●

DUAL-POINT CRYO-MARKER™

Cryogenic Waterproof Permanent Marker. High quality dual point permanent ink marker. Extra Fine (0.3mm) tip and fine (0.7mm) tip- Two pens in 1. Durable, heavy duty long lasting nibs (withstand rough treatment and provide delicate touch on various materials). Excellent performance on plastic, glass, metal, cloth and other surfaces. Can write on PTFE surfaces. Compatible with Lab-TAG™ labels. Sold per each.

Available Colors: Black, Blue, Red, Green, Purple, Orange, Brown and Yellow.

22315-2 through 22315-9 Dual-Point CRYO-MARKER™ each
 22315-10 Dual-Point CRYO-MARKER™, Assorted Colors.....pkg/8



www.tedpella.com sales@tedpella.com 800-237-3526

LABORATORY MARKING PENS

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS
www.tedpella.com

PERMANENT LABORATORY AND CRYO MARKING PENS



MICRO-TIP CRYO-MARKER™

Available colors:

Permanent and waterproof (all colors are equally permanent and lasting as black color) Excellent performance on plastic, glass, metal, cloth and other surfaces. Easy to handle. For applications in bio-medical laboratories Compatible with Lab-TAG™ labels. Durable, heavy duty long lasting nibs (withstand rough treatment and provide delicate touch on various materials).

Available Colors: Black, Blue, Red

22315-2 through 22315-13 Micro-Tip Cryo-Marker™ each
 22315-14 Micro-Tip Cryo-Marker™pkg/3




CRYO-EM MARKER PEN

Available colors:

Particularly suitable for writing on cryo-EM grid boxes and cryo storage vials, these markers feel and write like a pen. There's no offensive smell, and the cap closes with a nice reassuring click. The ink dries fast and is waterproof. Available in four colors in packages of 5 or 10.

Available Colors: Black, Blue, Red, Green

161-5B through 161-5R Cryo-EM Marker Penpkg/5
 161-10B through 161-10R Cryo-EM Marker Penpkg/10



BIG-TIP WAXI CRYO-MARKER™

Available colors:

Solid Ink. Permanent, easy to use, opaque, fade and waterproof. Marks on most surfaces and dries in 5-7 minutes. Recap after use. Ink may be removed with alcohols. Twist the end of barrel to advance marking tip. Tip point will not retract. Thermal resistance of ink: from -196°C to +200°C (-320°F to +392°F). Sold per each.

Available Colors: Black, Blue, Red, Green, Purple, Orange, Brown, Yellow, Fluorescent Orange, Fluorescent Yellow, and Fluorescent Pink.

27315-15 through 22315-24 Big-Tip Waxi Cryo-Marker™ each
 22315-25 Big-Tip Waxi Cryo-Marker™ Assorted Colorspkg/10



NISSEN™ PEN

Available colors:

The Nissen™ Nylon Tip Paint Marker dispenses a fast-drying enamel paint.

Marks metal, plastic, glass, wood, paper; Withstands heat; Available in black and 4 lead-free colors.

27305-1 Nissen™ Pen, Black each
 27305-A Nissen™ Pens, Assorted Colors (4 Black, 2 Red, 2 Green, 2 Blue, 2 Orange).....box/12



www.tedpella.com sales@tedpella.com 800-237-3526