

SLIDE & CASSETTE STORAGE

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS
www.tedpella.com



NEW PELCO® Microslides Storage Boxes

Exceptional slide storage containers providing maximum protection for your microscope slides. Features index card in the lid with associated slot number for each slide and foam liner. Boxes are stackable and made from durable/impact resistant ABS. Hinged lid and unique sliding compression latch. Available in 3 sizes and 6 colors.

25 slides: 104 L x 94 W x 32mm H (4.10" x 3.70" x 1.26")

50 slides: 164 L x 110 W x 32mm H (6.46" x 4.33" x 1.26")

100 slides: 209 L x 172 W x 32mm H (8.23" x 6.77" x 1.26")

2306-1	PELCO® Microslides Storage Box for 25 Slides, ○ White	each
2306-2	PELCO® Microslides Storage Box for 25 Slides, ● Red	each
2306-3	PELCO® Microslides Storage Box for 25 Slides, ● Green	each
2306-4	PELCO® Microslides Storage Box for 25 Slides, ● Blue	each
2306-5	PELCO® Microslides Storage Box for 25 Slides, ● Yellow	each
2306-6	PELCO® Microslides Storage Box for 25 Slides, ● Black	each
2304-1	PELCO® Microslides Storage Box for 50 Slides, ○ White	each
2304-2	PELCO® Microslides Storage Box for 50 Slides, ● Red	each
2304-3	PELCO® Microslides Storage Box for 50 Slides, ● Green	each
2304-4	PELCO® Microslides Storage Box for 50 Slides, ● Blue	each
2304-5	PELCO® Microslides Storage Box for 50 Slides, ● Yellow	each
2304-6	PELCO® Microslides Storage Box for 50 Slides, ● Black	each
2302-1	PELCO® Microslides Storage Box for 100 Slides, ○ White	each
2302-2	PELCO® Microslides Storage Box for 100 Slides, ● Red	each
2302-3	PELCO® Microslides Storage Box for 100 Slides, ● Green	each
2302-4	PELCO® Microslides Storage Box for 100 Slides, ● Blue	each
2302-5	PELCO® Microslides Storage Box for 100 Slides, ● Yellow	each
2302-6	PELCO® Microslides Storage Box for 100 Slides, ● Black	each

TED PELLA, INC.
 Microscopy Products for Science and Industry

www.tedpella.com sales@tedpella.com 800-237-3526

SLIDE & CASSETTE STORAGE

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS
www.tedpella.com



CRYO-COMPATIBLE BOX

Unique and durable polycarbonate construction with stainless steel latch lock. Protects 3" x 1" (75 x 25mm) slides in ultra-low temperature ranges of -80°C to -196°C (-112°F to -321°F). 100-slide capacity with grooved slots to prevent slide-to-slide contact. Offers a permanent identification sheet in hinged lid, that corresponds to the numbering on the temperature-resistant foam lining in base. The boxes are stackable and can be autoclaved. Made of polycarbonate.

2106 Cryo Compatible Slide Box, Holds 100 Slides, Purple..... each



TALL SLIDE BOX

Made of sturdy ABS material, this box is ideal for storage and transport of tall 3 x 2" (75 x 50mm) microscope slides. The grooved slots hold the slides securely in place. The base contains a cork liner with location numbers for convenient inventory. Hinged lid with a secure latch lock closure. The lid interior contains a permanent inventory sheet for cataloging and identification.

2197 Tall Slide Box, Holds 100 Slides each



SLIDE-SETTE SAVER™

The Slide-Sette Saver™ offers a convenient single container for shipping a paraffin tissue block together with its corresponding slide. The sturdy molded polypropylene container nestles the paraffin tissue block securely while offering a safe residence for its companion slide. The Slide-Sette Saver™ provides a protective storage environment for both block and slide in one unit and safely eliminates lost or separated patient samples. Available in natural or black.

Also available: Slide-Sette Saver™ Storage Box. Offers an organized, long-term storage system for 20 Slide-Sette Saver™ units.

22524-NAT Slide-Sette Saver™, Natural..... pkg/25
 22524-BLK Slide-Sette Saver™, Black pkg/25
 22523 Slide-Sette Saver™ Storage Box, Natural..... each



ECONOMY SLIDE STORAGE UNIT

The economy slide storage unit holds up to 1425 standard 3" x 1" (75 x 25mm) microscope slides. The unit features a single drawer, with reinforced metal corners. Drawer front has a finished metal edge and metal-edged drawer pull, as well as a label area for organization and cataloging

contents. The internal partitions keep the slides neat and easy to retrieve. The sturdy heavy pressed cardboard construction has a white satin finish, measures 12-1/2" x 7" x 4-1/2" (31.8 x 17.8 x 11.4cm), and is an economical alternative to the bulky, heavy metal storage cabinets for the active and growing laboratory.

27156 Slide Storage Unit, Holds 1425 slides.....each



TISSUE-TEK® CASSETTE FILING CABINET

Provides permanent storage-reference unit for blocks embedded in cassettes or embedding rings. The Tissue-Tek® Cassette Filing Cabinet accommodates 1500 cassettes or 1000 embedding rings. One-piece construction eliminates hinges and doors for dust-free storage. Includes six impact-resistant plastic drawers with identification area. Sold in USA only.

Dimensions: 43.2 x 24.1 x 38.1cm (17" x 9.5" x 15")

27169 Tissue-Tek® Cassette Filing Cabinet each



ECONOMY CASSETTE STORAGE UNIT

The economy tissue cassette storage unit holds up to 360 standard histology cassettes. Each drawer contains metal reinforced corners, a finished metal edge front and metal-edged drawer pull to enhance support and stability.

Equipped with strong internal partitions for organization and effective cataloging of specimens. The sturdy, heavy pressed cardboard construction has a white satin finish providing a clean, dust-free environment in a convenient 12-1/2" x 7" x 4-1/2" design.

27155 Cassette Storage Unit, Holds 360 Cassettes.....each