

# **NEW MICROSLIDES SLIDE FOLDER RACK**

For Cardboard or Plastic Pop-Up 20 Slide Folders



**STACKABLE  
MODULAR  
DESIGN**

- Holds cardboard pop-up 20-slide folders
- Holds plastic pop-up 20-slide folders
- Modular design

## **MICROSLIDES SLIDE FOLDER RACK**

Newly engineered, Slide Folder Rack. Each tray can accommodate either cardboard pop-up 20-slide folders (2110 series) or plastic 20-slide folders (2113 series). Unique design directs the slide folder into the appropriate slot providing equal access to each folder. The Slide Folder Rack comes with five trays and a cover. Modular design allows additional trays to be added to the rack. 5 color-coded slide folder clasps that keep the flaps of the slide folders closed are included. Additional clasps are available (product number 2110-10).



2110-8	Slide Folder Rack with 5Trays and Cover.....	each
2110-9	Slide Folder RackTray, Package of two additionalTrays.....	each
2110-10	Clasps for Pop-up Slide folders, 5 Assorted Colors, 20 Clasps of each color per package.....	each

For additional information on Slide Folders and Folder Clasps please visit us at [www.tedpella.com/histology](http://www.tedpella.com/histology)

 **TED PELLA, INC.**  
Microscopy Products for Science and Industry

[www.tedpella.com](http://www.tedpella.com)   [sales@tedpella.com](mailto:sales@tedpella.com)   800-237-3526