

SPOONS & SPATULAS

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS
www.tedpella.com

MICRO SPOONS



Mirror stainless steel, spoon and flat spatula. Spoon 29mm L x 14.3mm W (1-1/8" L x 9/16" W). Flat 48mm L x 14.3mm W (1-7/8" L x 9/16" W). Overall length 229mm (9").

13500 Micro Spoon pkg/3



Stainless steel, narrow spoon. Spoon 25.4 mm L x 6.4mm W (1" L x 1/4" W). Plastic handle, overall length 165mm (6-1/2").

13501 Micro Spoon pkg/3



Gray PTFE coated stainless steel, spoon and flat spatula. Spoon 29mm L x 14.3mm (1-1/8" L x 9/16" W). Flat 48 L x 14.3mm W (1-7/8" L x 9/16" W). Overall length 229mm (9").

13521 Micro Spoon/Spatula, PTFE Coated each

MICRO SPATULAS



Stainless steel, narrow spoon and flat spatula. Spoon 12.7mm L x 3.2mm W (1/2" L x 1/8" W). Flat 31.8mm L x 4mm W (2-1/4 L x 5/32" W), overall length 154mm (6-1/16").

13504 Micro Spatula pkg/3



Gray PTFE coated flat rounded ends. 41mm L x 7mm W (1.61" L x 1/4" W). Overall length 195mm (7.68").

13523 Micro Spatula, PTFE Coated each



Square & radius flat ends, stainless steel mirror finish. Radius end 51mm L x 8mm W (2" L x 5/16" W). Flat end 52mm L x 8mm W (2-1/16" L x 5/16" W). Overall length 178mm (7").

13508 Micro Spatula pkg/3



Rounded and rounded tapered ends. Rounded end 45mm L x 7mm W (1-5/16" L x 1/4" W). Rounded tapered end 49mm L x 7mm W (1-15/16" L x 1/4" W). Stainless steel, 200mm L (7-7/8").

13510 Micro Spatula pkg/3

©Ted Pella, Inc. 10-28-2013, Printed in U.S.A.