

PLASTIC MEMBRANE BOXES

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS

www.tedpella.com

Store/Ship/Transfer Delicate, Sensitive Articles Between Soft Elastic Membranes

Polystyrene plastic boxes with micro thin, flexible polyurethane membranes top and bottom. The plastic membranes form around the enclosed object and hold it securely. Small products or large may be held firmly within various box sizes. Semiconductor, machine parts, electron microscopy parts or materials, samples, and many other applications.



139-56 Holding Silicon Chips



139-58 Holding Silicon Wafer



139-51



139-59



139-58



139-50



139-60



139-63



139-56



139-52



139-64



139-62

Directions for use.




1.

Open box by pulling back the latch with your thumb.



2.

With box open, place objects on the lower membrane.



3.

Close lid and membranes will come together, deforming around the objects. Lock closed by pushing latch forward.

See reverse side for complete ordering information.

TED PELLA, INC.
Microscopy Products for Science and Industry

www.tedpella.com sales@tedpella.com 800-237-3526

PLASTIC MEMBRANE BOXES

FOR PRODUCT DETAILS AND COMPLETE SELECTIONS
www.tedpella.com

Store/Ship/Transfer Delicate, Sensitive Articles Between Soft Elastic Membranes

Membrane Boxes Ordering Information

139-50	Two Part, Press Together, Outside Dimensions: 38 x 38 x 16mm H (1.52" x 1.52" x 0.67" H), Membrane Area: 31.8 x 31.8mm (1.25" x 1.25").....	pkg/12
139-51	Two Part, Two Clear Clips, Outside Dimensions: 82mm dia. x 40mm H (3.22" dia. x 1.57" H), Membrane Area: 69.9mm dia. (2.75" dia.)	each
139-52	Two Part, Two Clear Clips, Outside Dimensions: 82mm dia. x 17.5mm H (3.22" dia. x 0.69" H) Membrane Area: 69.9mm dia. (2.75" dia.)	each
139-59	Two Part, Press Together, Outside Dimensions: 71 x 71 x 12mm H (2.8" x 2.8" x 0.47" H) Membrane Area: 60 x 60mm (2.36" x 2.36").....	each
139-60	Two Part, Press Together, Outside Dimensions: 39.6 x 79.4 x 19mm H (1.56" x 3.125" x 0.75" H) Membrane Area: 25.4 x 66.7mm (1" x 2.625").....	each
139-61	Two Part, Snap Together/Press to Remove, Outside Dimensions: 50 x 160 x 30mm H (1.97" x 6.3" x 1.18" H) Membrane Area: 30 x 140mm (1.18" x 5.51").....	each
139-56	Double Hinged, One Yellow Clip, Outside Dimensions: 75 x 124 x 26.7mm H (3" x 4.87" x 1.05" H) Membrane Area: 60.5 x 114.3mm (2.38" x 4.5").....	each
139-62	Single Hinged, One Yellow Clip, Outside Dimensions, Assembled: 100 x 50 x 25mm H (4" x 2" x 1" H) Membrane Area: 30 x 90mm (1.18" x 3.54").....	each
139-63	Double Hinged, One Yellow Clip, Outside Dimensions, Assembled: 100 x 75 x 25mm H (4" x 3" x 1" H) Membrane Area: 60 x 90mm (2.36" x 3.54").....	each
139-64	Single Hinged, One Yellow Clip, Outside Dimensions, Assembled: 100 x 50 x 50mm H (4" x 2" x 2" H) Membrane Area: 30 x 90mm (1.18" x 3.54").....	each
139-66	Double Hinged, Two Yellow Clips, Outside Dimensions, Assembled: 197 x 197 x 49mm H (7.75" x 7.75" x 1.94" H) Membrane Area: 180 x 187mm (7" x 7.375").....	each
139-58	Double Hinged, One Yellow Clip, Outside Dimensions: 149 x 149 x 51.8mm H (5.87" x 5.87" x 2.04" H) Membrane Area: 133.4 x 138.2mm (5.25" x 5.44").....	each

Replacement Clips and Hinges

139-51-10	Clear Replacement Clips for Membrane Box 139-51 and 139-52	pkg/10
139-57	Yellow Replacement Clips for Hinged Boxes 139-56, 139-58, 139-63, 139-64, 139-66	pkg/10
139-54	Clear Hinges for Hinged Boxes 139-56, 139-58, 139-63, 139-64, 139-66.....	pkg/10